Technology for Vacuum Systems

Instructions for use

Exhaust waste vapor condenser Peltronic™

After sales service:

Contact your local dealer or call +49 9342 808-5500.

Trademark index:

VACUU•LAN® (US-Reg.No 3,704,401), VACUU•BUS™, VACUU•CONTROL™, chemistry-HYBRID™, Peltronic™, TURBO•MODE™, VARIO® (US-Reg.No 3,833,788), VARIO-SP™, VACUUBRAND® (US-Reg.No 3,733,388) and also the shown company logos are trademarks of VACUUBRAND GMBH + CO KG in Germany and/or other countries.

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EN

Attention: This manual is not available in all languages of the EU. The user must not operate the device if he does not understand this manual. In this case a technically correct translation of the complete manual has to be available. The manual must be completely read and understood before operation of the device and all required measures must be applied.  "Safety instructions for vacuum equipment"

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Contents

Safety information!	9
Important information!	9
General information	11
Intended use.....	11
Setting up and installing the equipment	12
Ambient conditions	14
Operating conditions	14
Safety during operation	16
Maintenance and repair.....	18
⚠ Important information: Equipment marking (ATEX)	20
Technical data	22
Exhaust waste vapor condenser Peltronic™	22
Wetted parts	22
Abbreviations.....	23
Condenser parts.....	23
Exhaust waste vapor condenser Peltronic™ with CVC 3000.....	25
Exhaust waste vapor condenser Peltronic™ with VNC 2 (E) (VARIO)	25
Exhaust waste vapor condenser Peltronic™ with level sensor	26
Use and operation	27
Assembling components	27
Mounting the exhaust waste vapor condenser Peltronic™ at VACUUBRAND pumping units	27
Prior to use	29
During operation.....	30
Notes concerning the use of the exhaust waste vapor condenser Peltronic™	30
Shutdown & storage.....	31
Accessories	33
Troubleshooting	34
Notes on return to the factory	35
Warranty	38
Health and safety clearance form	39
EC Declaration of Conformity of the Machinery	40

Safety information!

Important information!

WARNING

- ☞ Keep this manual complete and accessible to personnel at all times!
- ☞ Read this manual carefully before installing or operating the equipment. Observe the instructions contained in this manual.
- ☞ Do not modify the equipment without authorization.

NOTICE

This manual is an integral part of the equipment described therein. It describes the safe and proper use of the exhaust waste vapor condenser Peltronic™. Make operating personnel aware of dangers arising from the exhaust waste vapor condenser Peltronic™ and the processed substances. VACUUBRAND disclaims any liability for inappropriate use of the exhaust waste vapor condenser Peltronic™ and for damage from failure to follow instructions contained in this manual.

This manual is only to be used and distributed in its complete and original form. It is strictly the users' responsibility to check carefully the validity of this manual with respect to his product.

Manual-no.: 999166 / 02/08/2012

The following signal word panels and safety symbols are used throughout this manual:

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury and death.

➔ DANGER indicates a hazardous situation which, if not avoided, will result in death or serious injury.

☞ WARNING indicates a hazardous situation which, if not avoided, could result in death or serious injury.

• CAUTION indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

NOTICE is used to address practices not related to personal injury.

Caution! Hot surface!

Disconnect equipment from AC power.

Dispose of electronic equipment according to regulations.

Formatting used in this manual:

Note: The signal word panels in all sections of this manual always refer to all paragraphs of the same format (➔ / ☞ / • / plain text) following each signal word panel.

The document "Safety information for vacuum equipment" is part of this manual! Read the "Safety information for vacuum equipment" and observe the instructions contained therein!

General information

NOTICE

Remove all packing material from the packing box. Remove the product from its packing-box and retain all packaging until the equipment is inspected and tested. Remove the protective caps from the inlet and outlet ports and retain for future use. Inspect the equipment promptly and carefully.

If the equipment is damaged, notify the supplier and the carrier in writing within three days. Retain all packing material for inspection. State the item number of the product together with the order number and the supplier's invoice number. Failure to check and give notice of damage will void any and all warranty claims for those deficiencies.

Replace the protective caps if the equipment is not used immediately. Store the equipment in dry and non-corrosive conditions.

WARNING

☞ **Do not use any damaged equipment.**

Intended use

WARNING

☞ Do not use any system parts on humans or animals.

☞ Ensure that the individual components are only connected, combined and operated according to their design and as indicated in the instructions for use.

☞ Comply with all notes on correct vacuum and electrical connections; see section "Use and operation", pg. 27.

☞ The exhaust waste vapor condenser Peltronic™ is designed for **ambient temperatures** during operation between +50°F and +104°F (+10°C and +40°C). Periodically check maximum temperatures if installing the equipment in a cabinet or a housing. Make sure ventilation is adequate to maintain recommended operating temperature. Install an external automatic ventilation system if necessary.

If pumping hot process gases, make sure that the maximum permitted gas inlet temperature is not exceeded.

CAUTION

- Ensure that the exhaust waste vapor condenser Peltronic™, and any accessories in the flow path are chemically resistant to the processed substances prior to operation.

NOTICE

Use the equipment **only as intended**, that is, for condensation of emitted vapors at the outlet of vacuum systems and pumps (i.e. at atmospheric pressure). Any other use will automatically invalidate all warranty and liability claims. Remain aware of safety and risks.

Setting up and installing the equipment

DANGER

- ➔ Equipment must be connected only to a **suitable electrical supply** and a suitable ground point. As such, the plug must be plugged into an outlet that is properly grounded. Failure to connect the motor to ground may result in deadly electrical shock.

The supply cable may be fitted with a molded European IEC plug or a plug suitable for your local electrical supply. The cable contains wires color coded as follows: green or green and yellow: ground; blue or white: neutral; brown or black: hot.

- ☞ Do not permit any **uncontrolled pressurizing**. Make sure that the exhaust pipeline cannot become blocked. If there is an exhaust isolation valve, make sure that you cannot operate the equipment with the valve closed to avoid a **risk of bursting!**
- ☞ Always provide a free and pressureless exhaust outlet to avoid risk of bursting.
- ☞ Keep the electrical power cord away from heated surfaces.

CAUTION

- Provide a firm, level platform for the equipment. Ensure a stable position of the exhaust waste vapor condenser Peltronic™. Check that all fittings are secure.
- Comply with **maximum permissible pressures** at inlet and outlet and with **maximum permissible pressure differences** between inlet and outlet. See section "Technical data", pg. 22. Do not operate the exhaust waste vapor condenser Peltronic™ with overpressure at the inlet.
- Check the power source and the condenser's rating plate to be sure that the power source and the equipment match in voltage, phase, and frequency.
- Make sure the fan's air supply is adequate, especially if the exhaust waste vapor condenser Peltronic™ is installed in an enclosure. Do not block ventilation slots.
- Connect hoses gas tight at inlet and outlet of the exhaust waste vapor condenser Peltronic™.
- **Note:** Flexible elements will shrink when evacuated.

NOTICE

Keep a minimum distance of 8 in (20 cm) between the cooling fan and surrounding items (e.g., housing, walls, etc.). Check ventilation slots regularly for dust/dirt. Clean ventilation slots if necessary to avoid a reduction of ventilation.

Allow the equipment to equilibrate to ambient temperature if you bring it from cold environment into a room prior to operation. Notice if there is water condensation on cold surfaces.

Use only hoses at the inlet and outlet of the condenser with an inner diameter at least as large as the diameter of the condenser's connections (to avoid overpressure, and reduction of throughput).

Comply with all **applicable and relevant safety requirements** (regulations and guidelines). **Implement the required actions and adopt suitable safety measures.**

Ambient conditions

DANGER

➔ Do not reach for this product if it has fallen into liquid. There is a risk of deadly electrical shock. Unplug the system immediately.

WARNING

☞ Do not use this product in an area where it can fall or be pulled into water or other liquids.

☞ Pay attention to the **permissible maximum ambient temperatures** (see "Technical data", pg. 22) and **gas temperatures**.

CAUTION

- Adopt suitable measures in case of differences from recommended conditions, e.g., using the equipment outdoors, installation in altitudes of more than 3300 ft (1000 m) above mean sea level, conductive pollution or external condensation on the exhaust waste vapor condenser Peltronic™.
- Do not operate this product near flames.

NOTICE

To the best of our knowledge the equipment is in compliance with the requirements of the applicable EC-directives and harmonized standards (see "Declaration of Conformity") with regard to design, type and model. Directive EN 61010-1 gives in detail the conditions under which the equipment can be operated safely (see also IP degree of protection, "Technical data", pg. 22).

Operating conditions

DANGER

➔ The exhaust waste vapor condenser Peltronic™ is not approved for operation in potentially explosive atmospheres. **Do not operate the exhaust waste vapor**

condenser Peltronic™ in potentially explosive atmospheres.

- ➔ Equipment **without the "Ex" mark on the rating plate is not approved for the condensation of potentially explosive atmospheres. Do not condense potentially explosive atmospheres with this equipment.**
- ➔ Equipment **bearing the "Ex" mark on its rating plate is approved for the condensation of potentially explosive atmospheres** according to its classification **II 3G IIC T3 X** according to ATEX, but it is **not approved for operation in potentially explosive atmospheres** (see section "Ex Important information: Equipment marking (ATEX)", pg. 20).

CAUTION

- **Do not condense** substances which may form **deposits** inside the exhaust waste vapor condenser Peltronic™. The exhaust waste vapor condenser Peltronic™ is not suitable for condensing substances which may form deposits inside the condenser. Deposits and condensate in the condenser may lead to increased temperatures even to the point of exceeding the maximum permitted temperatures.
- Check the inlet and outlet of the exhaust waste vapor condenser Peltronic™, if there is a danger of forming **deposits** inside the condenser. Inspect the exhaust waste vapor condenser Peltronic™ regularly and clean if necessary. See section "Shutdown", pg. 31).
- **Consider interactions and chemical reactions of the condensing media.** Ensure that the materials of the condenser's wetted parts are compatible with the processed substances, see section "Technical data", pg. 22.
When changing the substances processed, we recommend purging the exhaust waste vapor condenser Peltronic™ with air or inert gas prior to changing the processed media. Purging the exhaust waste vapor

condenser Peltronic™ will clean it from residues and it will reduce the possibility of reactions of the processed substances with each other and with the condenser's materials.

Safety during operation

- ➔ Adopt suitable measures to prevent the release of dangerous, toxic, explosive, corrosive, noxious or polluting fluids, vapors and gases. To prevent any emission of such substances from the condenser outlet, install an appropriate collecting and disposal system and take protective action for exhaust waste vapor condenser Peltronic™ and environment.
- ➔ You must take suitable precautions to prevent any formation of explosive mixtures in the condenser or at its outlet. In case, hot surfaces or static electricity may ignite these mixtures.
- ➔ Drain appropriately or otherwise remove any potentially explosive mixtures at the outlet of the condenser, or dilute them with inert gas to non-explosive concentrations.
- ➔ Never operate this equipment if it has a damaged cord or plug.

- ☞ If the condenser is not working properly, has been dropped or has fallen into water, contact your service provider.
- ☞ Make sure that the exhaust pipeline cannot become blocked.
- ☞ Comply with applicable regulations when disposing of chemicals. Take into consideration that chemicals may be contaminated. Take adequate precautions to protect people from the effects of dangerous substances (chemicals, thermal decomposition products of fluoro-

elastomers). Use appropriate protective clothing and safety goggles.

- ☞ Use only **original manufacturer's spare parts and accessories**. Otherwise the safety and performance of the equipment, as well as the electromagnetic compatibility of the equipment might be reduced.
The CE mark may be voided if not using original manufacturer's spare parts.
- ☞ Interruption of the equipment (e.g., due to power failure), failure of connected components or of parts of the supply, or change in parameters must not be allowed to lead to dangerous conditions. In case of a leak in the manifold, processed substances might be released into the environment.
Comply with all notes regarding proper use of the equipment, as well as operation and maintenance guidance.
- ☞ The residual **leak rate of the equipment** might render possible an exchange of gas, albeit extremely slight, between the environment and the system.
Adopt suitable measures to prevent contamination of the processed substances or the environment.

CAUTION

- Pay attention to the safety symbol "hot surfaces" on the equipment. Hot parts may cause burns if touched. Adopt suitable measures to prevent any danger arising from hot surfaces. Ensure that hot surfaces of the condenser do not cause burns. Provide a suitable contact guard if necessary.

NOTICE

Prevent the backpressure of gases and the backflow of condensates.

Check the liquid level in the catchpot during operation. Check the liquid level regularly. Do not allow the catchpot to overflow. Drain catchpot in time to avoid overflow.
Provide appropriate protective measures to allow for the

possibility of failure and **malfunction**. The protective measures must also allow for the requirements of the respective application.

An **excess temperature protection** switches off the condenser in case of excess temperature.

Attention: The condenser will start automatically after sufficient cooling down. It is your responsibility to ensure that automatic start-up of the condenser will not lead to any dangerous condition. Provide appropriate safety measures.

Maintenance and repair

NOTICE

In order to comply with laws (occupational, health and safety regulations, safety at work law and regulations for environmental protection) equipment can only be returned when certain procedures (see section "**Notes on return to the factory**", pg. 35) are followed.

! DANGER

➔ **Never operate the equipment if parts are disassembled.**

➔ **Switch off the exhaust waste vapor condenser Peltronic™. Disconnect the electrical power cord and wait two minutes** before starting maintenance to allow the capacitors to discharge.

➔ **Note:** The exhaust waste vapor condenser Peltronic™ may be contaminated with chemicals, which have been processed during operation. Ensure that the exhaust waste vapor condenser Peltronic™ is completely decontaminated before maintenance commences.

! WARNING

☞ Take adequate precautions to protect people from the effects of dangerous substances if contamination has occurred. Use appropriate protective clothing, safety goggles and protective gloves.

- ☞ **Never operate a defective or damaged exhaust waste vapor condenser Peltronic™.**
- ☞ Isolate the exhaust waste vapor condenser Peltronic™ from the vacuum system before starting maintenance. Allow sufficient cooling of the condenser. Drain condensate, if applicable.

NOTICE

Ensure that **maintenance** is done only by suitably trained and supervised technicians. Ensure that the maintenance technician is familiar with the safety procedures, which relate to the products processed by the system.

 Important information: Equipment marking (ATEX)**VACUUBRAND equipment bearing mark (see rating plate)**

** II 3G IIC T3 X
Internal Atm. only
Tech. File Ref.: VAC-EX01**

and

VACUUBRAND equipment bearing mark (see rating plate)

 X see manual

For equipment labelled with „ X see manual“ the following classification according to Directive 94/9/EC (ATEX) is valid:  II 3G IIC T3 X, Internal Atm. only, Tech. File Ref.: VAC-EX01.

The classification II 3G IIC T3 X according to ATEX is only valid for the inner part (wetted part, gases or vapors to be condensed) of the equipment. The equipment is not suitable for use in external, potentially explosive atmospheres (environment).

The overall category of the equipment depends on the connected components. If the connected components do not comply with the classification of the VACUUBRAND equipment, the specified category of the VACUUBRAND equipment is no longer valid.

Equipment in category 3 is intended for connection at equipment in which during normal operation explosive atmospheres caused by gases, vapors or mists normally don't occur or, if they do occur, are likely to do so only infrequently and for a short period only.

Equipment in this category ensures the requisite level of protection during normal operation.

The operation of venting valves is only permitted if thereby explosive atmospheres normally don't occur in the interior of the equipment or, if they do occur, are likely to do so only infrequently and for a short period.

The equipment is marked with "X" (according to EN 13463-1), i. e. restrictions of the operation conditions:

- The equipment is designated for a low degree of mechanical stress and has to be installed in a way so that it cannot be damaged from outside.
The equipment has to be installed protected against shocks from the outside and against glass splinters in the event of breaking (implosion).
- The equipment is designated for an ambient temperature during operation of +10 to +40°C. The equipment is designated for a gas inlet temperature during normal operation of +10 to +195°C. If the equipment is connected to the outlet of a vacuum pump, the permissible maximum gas inlet temperature of the vacuum pump must be observed (40°C if processing potentially explosive atmospheres. If processing gases/vapors which are not potentially explosive, extended gas inlet temperatures apply, see instructions for use of the vacuum pump.).
Never exceed these ambient and gas inlet temperatures.

After any intervention at the equipment (e.g., repair / maintenance) the leak rate of the equipment has to be checked. The measured leak rate has to be low enough to ensure that no explosive atmospheres will occur in the interior of the equipment.

Attention: This manual is not available in all languages of the EU. The user must not operate the device if he does not understand this manual. In this case a technically correct translation of the complete manual has to be available. The manual must be completely read and understood before operation of the device and all required measures must be applied.

Technical data

Exhaust waste vapor condenser Peltronic™		
Cooling capacity at 70°F (21°C) ambient temperature	W	approx. 50
Control temperature	°F (°C)	50 (10)
Permissible ambient temperature at operation	°F (°C)	50 to 104 (+10 to +40)
Permissible range of supply voltage (including tolerance)		90-126 V~ / 180-254 V~ 50-60 Hz
Rated current at: 100-120 V~ 50-60 Hz 200-230 V~ 50-60 Hz	A A	1.6 0.7
Power draw (controlled)	W	7 - 160
Heat dissipation including condensation heat	W	0 - 200
Permissible range of inlet and outlet pressure (absolute)	psi (bar)	13 to 16 (0.9 to 1.1)
Degree of protection IEC 529		IP 20
Inlet		connection for PTFE tube I.D. 3/8", O.D. 5/16" (10 / 8 mm)
Outlet		connection for PTFE tube I.D. 3/8", O.D. 5/16" (10 / 8 mm) or hose nozzle for tubing I.D. 3/8" (hose nozzle DN 10 mm)
Dimensions L x W x H approx.	in (mm)	6.9 x 7.0 x 15.4 (175 x 179 x 392)
Weight	lbs. (kg)	9.5 (4.3)
Volume of catchpot	quarts (ml)	0.52 (500)

Wetted parts

Components	Wetted materials
Fittings, hose nozzle	ETFE / ECTFE
Cooling surface	PP / PFA
Catchpot	Borosilicate glass (coated)

We reserve the right for technical modification without prior notice!


Abbreviations

ETFE:	Ethylene/Tetrafluoroethylene
ECTFE:	Ethylene/Chlorotrifluoroethylene
PFA:	Perfluoroalkoxylalkane
PP:	Polypropylene

Condenser parts

Position	Component
1	Inlet
2	Outlet
3	Mains connection
4	ON/OFF switch
5	LED
6	Connection plug of the VACUU • BUS line to the CVC 3000 or VNC 2 (VARIO) controller
7	Jack for connection of VACUU•BUS component (e.g., level sensor or valve or VMS module)
8	Catchpot

Exhaust waste vapor condenser Peltronic™


Exhaust waste vapor condenser Peltronic™ with CVC 3000


When switching on the controller, a connected exhaust waste vapor condenser Peltronic™ is detected automatically and supervised until the controller is switched off. The condenser is displayed by a symbol on the LCD of the controller with the actual condenser temperature shown in °C. The exhaust waste vapor condenser Peltronic™ is switched off simultaneously with the controller.

In case of fault or if removing the connection cable while the controller is switched on, a flashing warning triangle and a flashing condenser symbol (without temperature reading) are displayed. If "Sound" is set to "On", the controller sounds eleven blips.

Display of CVC 3000 during normal operation

condenser symbol
with condenser
temperature in °C

Display of CVC 3000
in case of fault

Exhaust waste vapor condenser Peltronic™ with VNC 2 (E) (VARIO)

When switching on the controller, a connected exhaust waste vapor condenser Peltronic™ is detected automatically and supervised until the controller is switched off.

The exhaust waste vapor condenser Peltronic™ is switched off simultaneously with the controller.

In case of fault or if removing the connection cable while the controller is switched on, "Ext. Err." ("external error") is displayed. If "Sound" is set to "On", the controller sounds additionally one blip.

Display of VNC 2 in case of fault

1013mbar
Ext. Err.

Exhaust waste vapor condenser Peltronic™ with level sensor

In case of high ambient humidity condensation may possibly occur on the cold surfaces of the exhaust waste vapor condenser Peltronic™. If condensation water is dripping onto the level sensor, an accidental level sensor alarm might be triggered. When installing a level sensor at the exhaust waste vapor condenser Peltronic™, pay attention to install it in a way that no potential condensation water can drip onto the level sensor. Avoid humidity on the sensor surface.


Valves and pressure transducers connected to the exhaust waste vapor condenser Peltronic™ are handled just as well as components connected directly to the controller.

NOTICE

For optimal electromagnetic compatibility assemble a ferrite at all VACUU•BUS connection cables (cat. no.: 612571).

Use and operation

Assembling components

Catchpot

- ➡ Assemble the catchpot at the condenser using the joint clip.

Separation of condensate:

The round bottom flask is coated with a protective layer to prevent disintegration in case of breakage.


- ☞ Separation of condensate: Remove joint clip, remove catchpot and drain condensate.

Important: Comply with applicable regulations when disposing solvents. Reuse if possible, purify if contaminated.

Mounting the exhaust waste vapor condenser Peltronic™ at VACUUBRAND pumping units

(not for VACUUBRAND NT pumping units with compact exhaust waste vapor condenser)

- ☞ For mounting the exhaust waste vapor condenser Peltronic™ at **VACUUBRAND NT pumping units with compact exhaust waste vapor condenser**, an assembly kit (order-no. 699935) is required.


Mechanical fastening:

- ➡ Dismount water cooled exhaust waste vapor condenser.
- ➡ Screw the exhaust waste vapor condenser Peltronic™ together with the holding block with two socket-head screws to the holder at the pump support.

Connecting the outlet of the pump to the exhaust waste vapor condenser Peltronic™:

- ☞ Establish connection with the formed PTFE hose included in delivery.
- ➡ Slip union nuts on hose. Fasten the PTFE hose with union nuts at the


tube fittings. If necessary widen the hose ends a bit by heating (538°F (280°C) at maximum).

☞ Depending on the pumping unit the formed PTFE hose has to be cut accordingly.

Attention: Do not crush the PTFE hose, in case use hose cutter (order number: 637257).

☞ Pumping units comprising a diaphragm pump with four pump heads and PC 2002: Assemble formed hose as delivered (1).

☞ Pumping units comprising diaphragm pumps with two pump heads (except PC 2002): Cut angled hose 90 mm x 60 mm (2).

☞ MD 1C AK + EK: Cut straight piece of hose. Length of 91 mm.

☞ PC 2001: Cut straight piece of hose. Length: 85 mm. Mechanical fastening only possible with special holding plate. Operation is also possible with free-standing condenser.

Hose cutter 637257

⚠ **Attention: The blade of the hose cutter is very sharp. Risk of injury!**

Prior to use

WARNING

- ☞ **Maximum ambient temperature:** 104 °F (40 °C)
Attention: High ambient temperature reduces the cooling power of the condenser. Take into consideration heating-up of the ambient air due to the vacuum pump and the exhaust waste vapor condenser Peltronic™.
- ☞ Check the maximum temperatures, if installing the condenser in a cabinet or a housing, or if the ambient temperature is elevated. Make sure ventilation is adequate.
- ☞ **Assemble exhaust waste vapor condenser Peltronic™ only at the outlet of a vacuum system (at atmospheric pressure).**

CAUTION

- Reduce the transmission of vibration. Prevent mechanical load due to rigid pipelines. Insert elastic hoses or flexible elements as couplings.
Note: Flexible elements will compress or flatten when evacuated if not designed for use under vacuum.
- If the equipment is installed at an altitude of more than 3300 ft (1000 m) above mean sea level, check compatibility with applicable safety requirements, especially IEC 60034. There is a risk of overheating due to insufficient cooling.
- If a power failure constitutes a potential source of danger, take appropriate safety measures.

NOTICE

Use connecting hoses with large diameter and keep them as short as possible to avoid flow losses. When assembling, ensure **tightness**. After assembly, check the whole system for leaks. Keep a distance of minimum 8 in (20 cm) between fans and adjacent equipment or casework.

Use the equipment **only as intended**, that is, for condensation of emitted vapors at the outlet of vacuum systems and pumps (i.e. at atmospheric pressure).

During operation

! DANGER

- ➔ Adopt suitable measures to prevent the formation of explosive or flammable mixtures, use inert gas for venting if necessary.
- ➔ Adopt suitable measures to prevent the release of dangerous, toxic, explosive, corrosive, noxious or polluting fluids, vapors and gases when disposing of condensates.


! WARNING

- ☞ **Maximum permissible operation pressure: 13 - 16 psi (0.9 - 1.1 bar) (absolute).**
Never block the gas outlet.

! CAUTION

- Make sure the fan's air supply is adequate. Clean ventilation slots, if necessary (e.g., remove dust or clean by suction).
- Prior to every use: Inspect the condenser for faults. There must be no damage. Do not use damaged components.

Notes concerning the use of the exhaust waste vapor condenser Peltronic™

The **exhaust waste vapor condenser Peltronic™** enables an efficient condensation of the pumped vapors at the outlet. The condenser's temperature is automatically controlled to +50°F (+10°C).

☞ Next to 100% solvent recovery.

- ☞ Never block the gas outlet. The exhaust hose has always to be unobstructed and without back pressure to enable an unhindered discharge of gases.
- ☞ Adopt suitable measures to prevent the release of dangerous or polluting fluids. To prevent any emission of such substances from the condenser outlet, install an appropriate collecting and disposal system.

NOTICE

The electronic exhaust waste vapor condenser Peltronic™ features a LED indicating the condenser's operation status (see also section "Troubleshooting").

Continuous green light: Condenser at operating temperature (cold side < +54°F (+12°C)).

Flashing orange light: Condenser temperature elevated (cold side > +54°F (+12°C)).

The condenser reaches its optimum cooling temperature even without load only about one minute after switching on. (LED: change from flashing orange to continuous green).

CAUTION

- **Continuous red light:** Excess temperature at the condenser's hot side. The cooling of the exhaust waste vapor condenser Peltronic™ is switched off. The fan remains running at maximum speed.

As soon as the condenser's hot side has cooled down to a temperature below 122°F (50°C), the exhaust waste vapor condenser Peltronic™ restarts automatically. The LED will be continuously cycling red and orange respectively red and green (depending on the temperature of the cold side). The flashing red light indicates that a fault had occurred.

The fault indicator "red light" can only be reset by switching the condenser off/on. Determine and eliminate the cause of failure (e.g., gas inlet temperature too high, gas load too high, ambient temperature elevated).

Shutdown & storage

NOTICE**Short-term:**

Remove catchpot and drain condensate.

Important: Comply with regulations when disposing of solvents/condensates. Recycle if possible; purify if contaminated.

Has the exhaust waste vapor condenser Peltronic™ been

exposed to media which may damage the condenser's materials or form **deposits**?

Check the exhaust waste vapor condenser Peltronic™ (check inlet and outlet) and clean, if necessary.

Clean a contaminated exhaust waste vapor condenser Peltronic™ as follows:

- Fill the condenser with a solvent (e.g., benzene) and allow sufficient cleaning time. Observe all regulations concerning usage and disposal of solvents!
- Drain the solvent and dispose of in accordance with regulations. Repeat cleaning if necessary.
- Rinse the condenser several times with alcohol in order to remove all solvent residues.
- Allow the emission condenser to dry.

Long-term:

- Take measures as described above regarding short-term shutdown.
- Separate the exhaust waste vapor condenser Peltronic™ from the application.
- Close inlet and outlet ports (e.g., with transport caps).
- Store the exhaust waste vapor condenser Peltronic™ under dry conditions.

Accessories

Assembly kit exhaust waste vapor condenser Peltronic™	699935
at VACUUBRAND NT pumping units with compact exhaust waste vapor condenser	
Straight fitting for PTFE tube DN 10/8 mm (at inlet or outlet)	638778 + 639796
PTFE tubing DN 10/8 mm, sold by meter	638644
Vacuum controller CVC 3000	683160
Vacuum controller VNC 2	683070
Vacuum controller VNC 2 VARIO	683080
Vacuum controller VNC 2 E	683090
Level sensor	699908
Pressure transducer VSK 3000, capacitive, ceramic diaphragm sensor 1080-0.1 mbar.....	636657
In-line valve VV-B 6, 24 V=	674290
In-line valve VV-B 6C, 24 V=	674291
In-line valve VV-B 15C, KF 16, 24 V=	674210
In-line valve VV-B 15C, KF 25, 24 V=	674215
Coolant valve VKW-B, 24 V=	674220
Venting valve VBM-B / KF 16, 24 V=	674217
Y-type adapter VACUU•BUS	636656
Extension cable VACUU•BUS, 2m	612552
VMS adapter cable (to connect Vacuum-Management-System VMS Module A)	636655
VMS-Module A, 100-230 V, 3,5 A, CEE	676001
VMS-Module A, 100-230 V, 3,5 A, UK	676002
VMS-Module A, 100-230 V, 3,5 A, CH	676003
VMS-Module A, 100-230 V, 3,5 A, US	676004
VMS-Module A, 100-230 V, 8 A, CEE	676006
VMS-Module A, 100-230 V, 8 A, UK	676007
VMS-Module A, 100-230 V, 8 A, CH	676008
VMS-Module A, 100-230 V, 8 A, US	676009
Spare parts:	
Catchpot 0.52 qt (500 ml), coated	638497

Troubleshooting

Fault	Possible cause	Remedy
<input type="checkbox"/> LED at the condenser is off, no fan noise.	<ul style="list-style-type: none"> ➔ Exhaust waste vapor condenser Peltronic™ not switched on or controller switched off? ➔ Electrical power cord not plugged in, electrical supply failure? ➔ Exhaust waste vapor condenser Peltronic™ defective? 	<ul style="list-style-type: none"> ✓ Switch on condenser and/or controller. ✓ Plug in power cord. Check fuse. ✓ Contact local distributor.
<input type="checkbox"/> LED is flashing orange, no change to continuous green light.	<ul style="list-style-type: none"> ➔ Very high condensation power? ➔ Ambient temperature too high? 	<ul style="list-style-type: none"> ✓ Check process parameters (e.g., reduce pumping speed of the vacuum pump, increase cooling at the vacuum side of the system). ✓ Comply with maximum permissible ambient temperature.
<input type="checkbox"/> Continuous red light at the condenser.	<ul style="list-style-type: none"> ➔ Ambient temperature much too high? ➔ Very high condensation power and ambient temperature too high? ➔ Exhaust waste vapor condenser Peltronic™ defective? 	<ul style="list-style-type: none"> ✓ Comply with maximum permissible ambient temperature. ✓ Check process parameters (e.g., reduce pumping speed of the vacuum pump, increase cooling at the vacuum side of the system). Comply with maximum permissible ambient temperature. ✓ Contact local distributor.

Notes on return to the factory

Repair - return - DAkkS calibration

NOTICE

Safety and health of our staff, laws and regulations regarding the handling of dangerous goods, occupational health and safety regulations and regulations regarding safe disposal of waste require that for all pumps and other products, the “**Health and safety clearance form**“, see pg. 117, must be sent to our office fully completed and signed before any equipment is shipped to the authorized service center.

Fax or mail a completed copy of the health and safety clearance form to us in advance. The declaration must arrive before the equipment. **Enclose a second completed copy with the product.** If the equipment is contaminated, you must notify the carrier.

No repair / DAkkS calibration is possible unless the correctly completed form is returned. Inevitably, there will be a delay in processing the equipment if information is missing, or if this procedure is not followed.

CAUTION

If the product has come in contact with chemicals, radioactive substances or other substances dangerous to health or environment, the product must be decontaminated **prior to sending it back to the service center.**

- Return the product to us disassembled and cleaned and accompanied by a certificate verifying decontamination or
- Contact an industrial cleaning and **decontamination service** directly or
- Authorize us to send the product to an industrial cleaning facility **at your expense.**

To expedite repair and to reduce costs, please enclose a detailed description of the problem and the product's operating conditions with every product returned for repair

We submit **repair quotations** only on request and always at the customer's expense. If an order is placed, the costs incurred for problem diagnosis are offset from the costs for repair or from the purchase price, if the customer prefers to buy a new product instead of repairing the defective one.

- **If you do not wish a repair on the basis of our quotation, the equipment may be returned to you disassembled and at your expense.**

In many cases, the **components must be cleaned in the factory** prior to repair.

For cleaning we use an environmentally friendly water-based process. Unfortunately the combined attack of elevated temperature, cleaning agent, ultrasonic treatment and mechanical stress (from pressurized water) may result in damage to the paint. Please mark in the health and safety clearance form, if you wish a **repaint at your expense** just in case such a damage should occur.

We will also replace parts for cosmetic reasons at your request and at your expense.

NOTICE

Before returning the equipment, ensure that (if applicable):

- Oil sealed pumps: Oil has been drained and an adequate quantity of fresh oil has been filled in to protect against corrosion. Dispose according to regulations.
- Equipment has been cleaned and/or decontaminated (inside and outside).
- All inlet and outlet ports have been capped.
- Equipment has been properly packed, (if necessary, please order original packaging materials at your cost), marked appropriately and the carrier has been notified of any possible contamination.
- The completed health and safety clearance form is enclosed.

We thank you in advance for your understanding of the necessity for these measures that protect our employees, and ensure that your pump is protected in shipment.

Scrapping and waste disposal:

Dispose of the equipment and any components removed from it safely in accordance with all local and national safety and environmental requirements. Particular care must be taken with components and waste oil which have been contaminated with dangerous substances from your processes. Do not incinerate fluoroelastomer seals and O-rings.

- You may authorize us to dispose of the equipment **at your expense.**

Warranty

VACUUBRAND shall be liable for insuring that this product, including any agreed installation, has been free of defects at the time of the transfer of risk.

VACUUBRAND shall not be liable for the consequences of improper handling, use, servicing or operation of this product or the consequences of normal wear and tear of wearing parts such as diaphragms, seals, valves, vanes, condensers, oil and the breakage of glass or ceramic parts, for the consequences of chemical, electrochemical or electrical influences or the failure to follow the instructions in this manual.

Claims for defects against VACUUBRAND shall be limited to one year from delivery. The same shall apply to claims for damages irrespective of legal grounds.

For further information on general terms and conditions refer to www.vacuubrand.com.

Health and safety clearance form

Devices will not be accepted for any handling before we have received this declaration.
Please read and comply with "Notes on return to the factory".
Oil filled pumps: Drain oil prior to shipping absolutely!

1. Device (Model): 2. Serial no.:

3. Reason for return / malfunction:

.....

4. Has the device been used in a copper process step (e.g., semiconductor production).

yes no

5. Substances (gases, liquids, solids) in contact with the device / which have been pumped:

.....

.....

.....

6. Prior to return to the factory the device has been decontaminated.

yes no

Description of the decontamination method and the test / verification procedure:

.....

7. The device is free of hazardous, harmful substances.

yes no

8. Protective measures required for VACUUBRAND employees:

.....

9. If the paint is damaged, we wish a repaint or a replacement of parts for reason of appearance (repaint and replacement at customer's expense).

yes no

10. Legally binding declaration

We assure for the returned device that all substances, which have been in contact with the device are listed in section 5 and that the information is complete and that we have not withheld any information. We declare that all measures - where applicable - have been taken listed in section "Return to the factory". By our signature below, we acknowledge that we accept liability for any damage caused by providing incomplete or incorrect information and that we shall indemnify VACUUBRAND from any claims as regards damages from third parties. We are aware that as expressed in § 823 BGB (Public Law Code of Germany) we are directly liable for injuries or damages suffered by third parties, particularly VACUUBRAND employees occupied with handling/repairing the product.
Shipping of the device must take place according to regulations.

Name: Signature:

Job title: Company's seal:

Date:

Release for repair grant by VACUUBRAND (date / signature):

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Hereby the manufacturer declares that the machinery is in conformity with the directive 2006/95/EC.

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Disclaimer: Our technical literature is only intended to inform our customer. The applicability of general empirical values and results obtained under lab conditions to your specific operations depends on a number of factors beyond our control. It is, therefore, strictly the users' responsibility to very carefully check the application of these data to their specific requirements. No claims arising from the information provided in this literature will, consequently, be entertained.

Alfred-Zippe-Str. 4 · 97877 Wertheim / Germany
T +49 9342 808-0 · F +49 9342 808-5555
info@vacuubrand.com · www.vacuubrand.com

VACUUBRAND GMBH + CO KG
- Technology for Vacuum Systems -
© 2012 VACUUBRAND GMBH + CO KG Printed in Germany